



EstheticoTM

PRODUCT CATALOG



Esthetico®

ABOUT US



OUR VALUES

Esthetico constantly makes an effort to understand consumers to provide them with utmost satisfaction.

Esthetico team continuously monitors consumer trends and reflect them into it's products.



DEFINITION

Esthetico logo symbolizes internal and external confidence in looks and all walks of life.



SPIRIT

Expertise in management based on respect & value creation for customers to realize a safer beauty life.



IDENTITY

Esthetico is a global leader in health and beauty devices that is focused on achieving an utmost customer's satisfaction with the highest level of technology and quality in the medical device.



VISION

We will provide solutions and services and be a preferred global aesthetic company that brings smile in a sustainable way.



SLOGAN

You Trust, We Care



www.esthetico.co.kr



Esthetico®

CERTIFICATION



Venture for
Tomorrow



www.esthetico.co.kr

UCO-HIGH INTENSITY FOCUSED ULTRASOUND SYSTEM

Uco – System Introducing our High Intensity Focused Ultrasound (HIFU) System – the cutting-edge solution designed to address the diverse needs of achieving beautiful and flawless skin. This advanced system has been meticulously crafted to fulfil eight key conditions for optimal skin enhancement.

INDICATIONS:

HIFU Beauty Boost: Elevate Your Radiance!

Unleash the power of our High Intensity Focused Ultrasound (HIFU) System for a quick, non-surgical beauty makeover. Experience:

- **Face Lifting & Tightening** – Sculpt your features and defy gravity.
- **Wrinkle Reduction**– Smooth away lines and embrace a refined complexion.
- **Face Rejuvenation**– Awaken a youthful glow and rediscover your essence.
- **Eyebrow Lifting**– Lift your gaze and enhance your allure effortlessly.
- **Facial/Slim Type**
 1. **(HD-1.5mm)** Remove wrinkles around the eyes, eye brows, nose
 2. **(HD-3.0mm)** Remove nasolabial line and wrinkles in facial area, including forehead and cheeks
 3. **(HD-4.5mm)** Remove nasolabial line and wrinkles in cheeks, double chins, neck, chest
- **Body**
 1. **(HD-7mm)** Remove wrinkles around the eyes, eye brows, nose
 2. **(HD-10mm)** Sides, forearms, thighs, abdomen
 3. **(HD-13mm)** Abdomen, forearms, thighs, sides, hip

Technical Specification:

Energy Type	High Intensity Focused Ultrasound	
Power	AC 100-230 V 50/60 Hz	
Display	TFT 8 Inch Color LCD Touch Screen	
Energy Fluency	0.2~2J (0.1 step)	
Cartridge Type	Treatable Parts Of Body	
Facial/Slim Type	HD-1.5mm	Remove wrinkles around the eyes, eye brows, nose
	HD-3.0mm	Remove nasolabial line and wrinkles in facial area, including forehead and cheeks
	HD-4.5mm	Remove nasolabial line and wrinkles in cheeks, double chins, neck, chest
Body	HD-7mm	Remove wrinkles around the eyes, eye brows, nose
	HD-10mm	Sides, forearms, thighs, abdomen
	HD-13mm	Abdomen, forearms, thighs, sides, hip
Cartridge Shot	10,500 shot	
Spacing	1.0~2.0mm (0.5mm/step), 2.0~5.0mm(0.5mm/step)	
Length	5.0~28mm (1.0mm/step)	
Weight	5.7kg (excluding cart)	
Dimensions	370x310x850mm (WxDxH)	
Cart	Optional	
Cart Dimensions	350x360x850mm (WxDxH)	



MATRIX – CO2 FRACTIONAL LASER 20W

Matrix – Introducing our Matrix CO2 Fractional Laser 20W – an advanced solution for comprehensive skincare and aesthetic treatments. This cutting-edge device is designed to address a range of dermatological concerns with precision and efficiency.

INDICATIONS:

- **Scars remodeling** – Effectively targets and remodels scars for improved skin texture.
- **Skin renewing and resurfacing** – Promotes skin renewal, leaving a smoother and rejuvenated complexion
- **Wrinkle removal and skin tightening** – Reduces wrinkles and tightens the skin for a youthful appearance
- **Pigmentation removal** – Targets and eliminates unwanted pigmentation, restoring even skin tone
- **Skin excrescence** – Addresses skin excrescence with precision and accuracy
- **Acne treatment** – Provides a powerful solution for acne treatment, aiding in clearer skin
- **Vaginal treatment** – Offers specialized treatment for vaginal concerns, contributing to overall feminine wellness

*

Technical Specification:

Power	Max. 20W
Fractional Laser	Pulse Energy : 3 ~ 300 [mJ] Scan Area Size : upto 20 x 20 [mm]
Continuous Wave	Spot Size(30W) : 50mm Hand Piece (1.0 [mm]), 100mm Hand Piece (1.5 [mm]) Power : Max. 20W
Ultra Pulse	Frequency : 10 ~ 1,000 [Hz] Duration : 100 ~ 2,000 [us]
Display	12.1" LCD
Aiming Beam	Diode Laser Beam 635nm
Laser Delivery	7Joint Articulated Arm
Input Voltage	200 ~ 240 [VAC], 50/60 [Hz]
Weight	40 [kg]
Dimension	460 x 430 x 1,170 [mm] (W x D x H)





GRAND – CO2 FRACTIONAL LASER 30W

Grand – Lets dive into the journey of cosmic realm of each indication for our Grand CO2 Fractional Laser System:

Radiate Cosmic Beauty: CO2 Laser Magic!

Experience a stellar makeover with our Grand CO2 Fractional Laser System. Shrink pores, smooth wrinkles, and transform scars into cosmic wonders. Find comfort with hemorrhoid relief, clarity with condyloma eradication, and sparkle with spot and lentigines dissolution. Watch corns vanish into the cosmic abyss. Elevate intimacy with vaginal tightening and whitening. Unveil a universe of radiant possibilities with each pulse. Let Beauty be redefined.

INDICATIONS:

Large Pore Refinement

- Pore Nebula Resculpt.
- Witness a celestial makeover as our laser system sculpts pores into miniature galaxies, creating a silky-smooth cosmic

Wrinkle Smoothing

- Timeless Celestial Elixir.
- Immerse yourself in the enchanting elixir of our laser's touch, as it orchestrates a cosmic ballet, smoothing away wrinkles and revealing the celestial glow within.

Acne Scar Transformation

- Scar Constellation Dissolution.
- Watch scars transform into a mesmerising celestial dance, fading away like distant constellations in the night sky.

Hemorrhoid Relief

- Nebular Comfort Waves
- Experience the soothing waves of comfort as our laser system creates nebular relief for hemorrhoidal concerns.

Condyloma Clarity

- Galactic Eradication Ripples
- Illuminate your skin with ripples of clarity as unwanted condylomas vanish, leaving behind a radiant complexion.

Spot Illumination

- Cosmic Sparkle Revelation
- Behold the cosmic sparkle as spots dissolve into radiant luminosity, unveiling a dazzling, even-toned skin.





Lentigines Radiance

- Stellar Lentigo Luster
- Shine with the luminous glow of the stars as lentigines surrender to the laser's transformative radiance.

Corn Dissipation

- Nebula Corn Vanishing Act
- Marvel at the cosmic vanishing act as corns dissipate into the galactic abyss, restoring a sense of celestial comfort.

Vaginal Tightening

- Celestial Intimacy Renewal
- Embark on an intimate journey of renewal as our laser system orchestrates a cosmic dance, tightening and revitalising.

Luminescence Whitening

- Venusian Glow Ritual
- Illuminate your inner beauty with a celestial ritual, imparting a whitening glow that transcends earthly radiance.

Technical Specification

Power	Max. 30W
Fractional Laser	Pulse Energy : 3 ~ 300 [mJ] Scan Area Size : 3 x 3 ~ 20 x 20 [mm] Scan Area Shape : Square, Rectangular
Continues Wave	Spot Size(30W) : 50mm Hand Piece (1.0 [mm]), 100mm Hand Piece (1.5 [mm]) Power : Max. 30W
Ultra Pulse	Frequency : 10 ~ 1,000 [Hz] Duration : 100 ~ 2,000 [us]
Display	8" LCD
Aiming Beam	Diode Laser Beam 650nm (Light Adjustable)
Laser Delivery	7Joint Articulated Arm
Input Voltage	200 ~ 240 [VAC], 50/60 [Hz]
Weight	62 [kg]
Dimension	450 x 510 x 1,150 [mm] (W x D x H)



ACTIV-Q-SWITCH ND YAG LASER SYSTEM 1500MJ

Activ – Embark on a luminous odyssey with our Activ Q-Switch ND Yag Laser System 1500mj – the epitome of skincare elegance! Let the cosmos of wavelengths, from the standard 1064 nm and 532nm to optional frequencies like 585 nm, 650nm, 755nm, and Fractional Handpiece versatility, redefine your skincare experience.

INDICATIONS:

- **Skin Rejuvenation** – Uncover a radiant complexion, bathing in the glow of rejuvenated skin.
- **Freckle Freedom** – Bid farewell to freckles as they gracefully fade into the canvas of flawless skin.
- **Brown Spots**– Watch brown spots dissolve, leaving behind a seamless complexion.
- **Seborrheic Keratosis**– Erase imperfections for a smooth, revitalized appearance.
- **Sun Spot Vanishing Act**– Witness sun spots disappear, revealing a brighter complexion.
- **Tattoo Erasure**– Basic color tattoos surrender to the precision of Q-Switch ND Yag Laser magic.
- **Nevus Harmony**– Tackle Fusco-Caeruleus, Zygomaticus, Ota, and Lto nevi with confidence.
- **Mongolian Spot Magic**– Illuminate your skin as Mongolian spots vanish into the cosmic glow.
- **Chloasma Serenity**– Achieve tranquility as chloasma dissolves, leaving a luminous complexion.
- **Pigmentation Perfection**– Unleash a canvas of flawless skin with the precision of our advanced laser system.
- **Carbon Peeling**– Reveal refreshed skin with the exfoliating touch of carbon peeling, adding another dimension to your skincare ritual.

Technical Specification:

Laser Type	Electro-optic nd yag Q-switched laser
Wavelength	(Standard) 1064nm, 532nm, (Optional) 585nm/ 650nm/ 755nm / Fractional Hand Piece (Optional)
Laser output method	Pulse electro-optic adjusting Q
Pulse Energy	1500mj (1064nm); 800mj (532nm)
Power	2000w
Width of pulse	6ns
Frequency	1-10Hz
Diameter of spots	2-10mm adjustable
Beam of light transmission	7 articular-arm of light guide transmitting
Indicator of aiming light	Red diode laser
Cooling manner	Closed-off water circulation+air
Net Weight	85kg
Voltage	AC200v+10% 50Hz, 110v+10% 60Hz





ACTIV PLUS – Q-SWITCH ND YAG LASER 2700MJ

Active Plus – Illuminate Your Skin’s Journey with Activ Plus – Q-Switch ND Yag Laser Elegance!

Introducing our Activ Plus-Q-Switch ND Yag Laser System, a beacon of advanced skincare. Immerse yourself in a world of wavelengths, including the standard 1064 nm and 532nm, with optional frequencies like 585 nm, 650nm, 755nm, and Fractional Handpiece versatility.

INDICATIONS:

- **Skin Rejuvenation** – Unveil a radiant complexion, basking in the glow of rejuvenated skin.
- **Freckle Freedom** – Bid farewell to freckles as they fade into the canvas of flawless skin.
- **Brown Spots**– Watch brown spots dissolve, leaving behind a seamless complexion.
- **Seborrheic Keratosis**– Erase imperfections for a smooth, revitalised appearance.
- **Sun Spot Vanishing Act**– Witness sun spots disappear, revealing a brighter complexion.
- **Tattoo Erasure**– Basic colour tattoos surrender to the precision of Q-Switch ND Yag Laser magic.
- **Nevus Harmony**– Tackle Fusco-Caeruleus, Zygomaticus, Ota, and Lto nevi with confidence.
- **Mongolian Spot Magic**– Illuminate your skin as Mongolian spots vanish into the cosmic glow.
- **Chloasma Serenity**– Achieve tranquillity as chloasma dissolves, leaving a luminous complexion.
- **Pigmentation Perfection**– Unleash a canvas of even skin tone, free from pigmentation concerns.
- **Carbon Peeling**– Reveal refreshed skin with the exfoliating touch of carbon peeling.

Technical Specification:

Wavelength	1064 [nm] / 532 [nm]
Operation Method	Q-switched Nd:YAG
Pulse Energy 1064 mode(Q-switched Mode)	Max. 2700 [mJ]
Pulse Energy 532 mode (Q-switched Mode)	Max. 900 [mJ]
Type of Pulse	Single Pulses : One whole pulse Double Pulses : Photoacoustic twin pulses Triple Pulses : Photoacoustic triple pulses
Pulse Width	5 ~ 10 [ns] (Q-switched Mode) / 500 [μs] (FR mode)
Spot Size	2 ~ 10 [mm]
Pulse Rate 1064, 532, FR Mode	1 ~ 10 [Hz]
Memory	5 User Programmable Memory
Optical Delivery	Articulated Arm
Aiming Beam	Diode 650 [nm] (Adjustable Brightness)
Display	10.4" LCD
Cooling	Closed Circuit Water to Air
Input Power	100 ~ 240 [VAC], 50/60 [Hz]
Dimension	315 x 825 x 1,022 [mm] (W x D x H)
Weight	90 [kg]



V1200-ULTRASONIC WAVE INNOVATION

V1200 – Incorporating ultrasonic waves, the machine delivers controlled energy to break down fat cells, allowing for a more comprehensive and efficient fat reduction process. This non-invasive method ensures precision in targeting localized fat deposits.

INDICATIONS:

- **Versatile Weight Loss**

Target stubborn fat, improves the skin texture and enhance overall body contour.

- **Non-Invasive Precision**

Achieve results without the need for surgery or invasive procedures.

- **Tailored Treatments**

Customised sessions to address specific body areas and individual weight loss goals.

- **Comfortable Sessions**

Enjoy a relaxing experience with controlled vacuum and pressure dynamics.

- **Visible Results**

Experience noticeable improvements in body shape, skin tone, and cellulite reduction.

Transform your weight loss journey with the Infrared RF Vacuum Roller Slimming Machine – a fusion of advanced technologies working harmoniously to sculpt your body and unveil the radiant and confident you.

*

Technical Specification:

Power of RF	1200w
Frequency of RF	2Mhz
Infrared laser Power	20w
Infrared laser wavelength	940nm
Ultrasound power	200W
Vacuum Indesity	10-90Kpa
Dimensions	470*500*1330mm
Weight	50KG





PULSE-PICO-SECOND LASER SYSTEM

Pulse – Enter the Quantum Realm of Flawless Beauty: Pico-Second Laser System
Unleash a revolution in skin care with our Pico-Second Laser System – a cosmic force of rejuvenation and transformation. Traverse the indications, each a celestial invitation to reveal your true radiance.

INDICATIONS:

- **Skin Rejuvenation** – Experience the revitalizing touch as our device unveils a radiant and youthful complexion.
- **Freckle** – Bid farewell to freckles as they gracefully fade away, leaving behind a flawless canvas.
- **Brown Spots** – Watch brown spots dissolve, revealing an even and luminous skin tone.
- **Seborrheic Keratosis** – Effortlessly erase imperfections for a smooth and revitalized appearance.
- **Sun Spots** – Witness the disappearance of sun spots, unveiling a brighter and more radiant complexion.
- **Tattoo Removal for Basic Colors** – Harness the power of precision for effective removal of basic color tattoos, achieving a clean slate.
- **Nevus** – Fusco-Caeruleus, Zygomaticus/Ota/Lto – Tackle various types of nevi with confidence, promoting harmonious skin.
- **Mongolian Spot** – Illuminate your skin as Mongolian spots vanish into the cosmic glow, leaving behind pure radiance.
- **Chloasma** – Achieve tranquility as chloasma dissolves, leaving a luminous and even skin tone.
- **Pigmentation** – Unleash a canvas of flawless skin by targeting and perfecting various pigmentation concerns.
- **Carbon Peeling** – Reveal refreshed and rejuvenated skin with the exfoliating touch of carbon.

Technical Specification:

Wavelength	Standard:1064nm/532nm. Optional:585nm/650nm/755nm
Energy	1500mj (1064nm) ; 800mj (532nm)
Peak Power	1064nm 1.33GW; 532nm 0.67GW
Frequency	1-10Hz
Zoom Spot Size	2-10mm Adjustable
Pulse Width	450ps
Beam Profile	Top Hat Beam
Light Guiding System	7 joints Arm
Aiming Beam	Diode 655 nm (Red)



ROMA-DIODE LASER 600/800W

Roma – Step into a new era of hair removal excellence with our Roma Diode Laser 600/800W. This state-of-the-art device is a technical marvel, designed to redefine your expectations and elevate your aesthetic experience.

Elevate your aesthetic practice with Roma Diode Laser – where cutting-edge technology meets versatility, offering a multitude of features for optimal hair removal. Revolutionize your approach to skincare with a device that seamlessly combines innovation and effectiveness.

Technical Specification:

Light Source	Harness the power of a Diode Module for superior performance
Output Energy	Choose between 600 [W] with 6 bars or 800 [W] with 8 bars, offering versatility for various treatment needs
Shots	Enjoy an extensive lifespan with an impressive 10,000,000 shots, ensuring enduring efficacy
Wavelength	The 808 [nm] wavelength ensures precision targeting for effective hair removal
Fluence	Experience optimal energy delivery with 100 [J/cm ²] (600[W]) or 102.7 [J/cm ²] (800[W])
Pulse Duration	Customizable pulse duration ranging from 10 to 350 [ms] for tailored treatments
Spot Size (Sapphire Size)	Choose your spot size – 14 x 15 [mm] (600[W]) or 14.2 x 18.5 [mm] (800[W]) for precise application
Repetition Rate	Flexible repetition rates from 1 to 10 [Hz] for personalized treatment plans.
Input Voltage	Adaptable input voltage of 100-240 [VAC], 50/60 [Hz] for global compatibility
Dimension	Compact dimensions at 420 X 510 X 1,040 [mm] (W x D x H) for efficient use of space
Weight	A manageable weight of 50 [kg], ensuring portability without compromising on performance



RAZOR-DIODE LASER SYSTEM 500W

Razor – Unleash the future of hair removal with our state-of-the-art Razor Diode Laser System 500W, a revolutionary solution designed to redefine your aesthetic experience.

INDICATIONS:

- **Permanent Hair Removal Solution** – Bid farewell to unwanted hair with our advanced diode laser technology.
- **Easy to Operate** – User-friendly design ensures a seamless and straightforward operational experience.
- **Safe on All Skin Types** – Enjoy the confidence of safe and effective treatments across diverse skin tones.
- **Virtually Painless** – Experience comfort with virtually painless hair removal sessions.
- **Faster Delivery Speeds** – Swift and efficient, our laser system offers accelerated treatment delivery.
- **Efficient Full Body Treatments** – Achieve comprehensive results with efficient full-body treatments.
- **Ability to Treat Large or Small Areas** – Versatile capabilities for addressing various treatment areas.

- **Comfortable & Light Weight Hand Piece Design** – Ergonomic design for comfort during extended use.
- **Zero Consumables** – Eliminate the need for consumables, ensuring cost-effective and sustainable treatments.

Technical Specification:

Laser Type	Diode laser
Wavelength	3 in 1 Wavelength, or single 808nm or single 1064nm
Safety Level	Ensures a high standard of safety for every treatment
Pulse Width	Adjustable from 5-400ms
Frequency	Up to 10Hz (max)
Spot Size	12mm*10mm; (12mm*16mm optional)
Energy Density	Up to 120J/CM2 (max)
Refrigeration System	Built-in water cooling, air cooling, semiconductor cooling
Power	1300VA
Voltage	110V-240V AC, 50/60Hz
Cooling System	Triple cooling system for optimal performance
Net Weight	40kg
Dimension (WxDxH)	470mm*470*1100mm



EMS 5000- EMSCULPT – THE POWER OF MUSCLE SCULPTING

EMS 5000 – Unleash the future of body transformation with EMSCULPT, where muscle sculpting meets innovation. This cutting-edge technology harnesses non-invasive, highly focused electromagnetic waves, delving up to 7CM deep into muscle tissue. It's a symphony of precision, triggering motor neurons, and orchestrating super-strong contractions – an extreme workout at frequencies approaching 100% muscle fibre activity.

Technical Specification:

Power	5000W
RF Power	800w
Pulse width	310us±10%
SPulse shape	two-way wave
Screen	12.4inch tough screen
Frequency:	1-150hz
Voltage	AC220V±10% 10A50HZ,110V±10% 10A 60HZ
Cooling:	Air Cooling
intensity	7 Tesla
Handles	4 Handles





Esthetico®

THANK YOU

Contact Us

📍 Medical Device Complex Center, 200 Gieopdosi-ro, Jijong-myeon, Wonju-si, Gangwon-do, 26354, Republic of Korea

✉️ info@esthetico.co.kr

www.esthetico.co.kr